



Baxter

Floseal
HEMOSTATIC MATRIX

**UNMATCHED SPEED.
RELIABLE HEMOSTASIS.**



Please see Indication and Important Risk Information on page 12



PATIENT PROFILES ARE CHANGING, INCREASING THE RISK OF BLEEDING RELATED COMPLICATIONS

The surgical patient profile has changed over the years: increased age, comorbidities, and the increased utilization of antiplatelet/anticoagulant products may impact a patient's coagulation status.¹

Surgical bleeding is associated with an increased risk of morbidity and mortality across the spectrum of surgical specialties.² Bleeding-related complications have been found to occur in **nearly 30% of surgeries overall, with the rate approaching 50% for cardiac-related surgery.**²



2/3
of surgical patients are treated
with anticoagulants and/or
antiplatelets^{3*}

*Based on Premier Database



1 in 50
adults have liver disease, a
comorbidity that negatively
impacts coagulation status⁴



~30%
of the US population uses
low-dose aspirin⁵

Understand the difference to select the right product

Managing bleeding requires a patient-centric approach. Surgical teams benefit from an understanding of the different types of hemostatic agents available to supplement a surgeon's skills in achieving meticulous hemostasis.

A

ACTIVE HEMOSTATIC AGENTS¹

Function independently of the patient's ability to generate clotting factors to achieve hemostasis and are effective regardless of patient coagulation status in a wide range of bleeding.¹



Flowables
(thrombin + gelatin)



Fibrin sealant



Stand-alone
thrombin



Advanced patches with
PEG or fibrin sealants

P

PASSIVE HEMOSTATIC AGENTS¹

Aid in platelet aggregation and/or contact activation for platelet plug formation and rely on the patient's ability to generate clotting factors to convert fibrinogen to help form the fibrin clot.¹



Oxidized
cellulose



Collagens



Polysaccharide
spheres



Powders



Gelatin
sponges


IMPORTANCE OF CHOOSING THE **RIGHT** HEMOSTATIC AGENT

Not knowing which adjunctive hemostatic agent to select can have negative impacts^{1,6,7}

- Delays in hemostasis
- Excessive use of multiple hemostatic agents
- Increased intraoperative times
- Postoperative complications
- Extended hospital stay
- Increased costs

LEARN THE VIBE SCALE TOOL: SEVERITY OF BLEED INFORMS CHOICE OF ADJUNCTIVE HEMOSTATIC AGENT^{8,9}

The Validated Intraoperative Bleeding (VIBe) SCALE tool classifies bleeding severity by amount of blood loss, allowing surgical teams to effectively communicate the level of bleed.^{8,9} Passive hemostatic agents are appropriate for Grade 1 bleeds for patients with intact coagulation.¹ For patients with compromised coagulation, consider an active hemostatic agent that can address bleeds that are mild, moderate, or severe - regardless of the patient's coagulation status.¹




FLOSEAL MATRIX: FAST HEMOSTASIS AND PREPARATION^{10,12}

Hemostasis¹⁰

FLOSEAL Matrix has the potential to achieve hemostasis in as little as

60 seconds*

*For spinal/orthopedic specialty, median time to hemostasis was 1.5 min (95% CI, 1, 1.5)



Proprietary **Floseal** Matrix Smooth Granules¹¹

Preparation¹²

FLOSEAL Matrix can be prepared in approximately

1 minute


FLOSEAL MATRIX IS EFFECTIVE AND PROVEN^{10,14}

Effective¹⁰

FLOSEAL Matrix is effective across multiple surgical specialties at a rate of

96%*

*Hemostasis within 10 minutes - first lesion only



†Floseal Matrix vs. Gelfoam Hemostat + Thrombin. Median time to hemostasis 2 minutes vs. 6 minutes.

Proven^{10,14}

Proven hemostasis and relative reductions in blood transfusions in cardiac surgery of

54%†

†Transfusion of blood products for **Floseal** Matrix n=31 [28.2%] vs comparator group [**Surgical Nu-Knit** Hemostat or **Gelfoam** sponge n=63 [60.6%, p<0.001]].




SCAN TO WATCH THE MECHANISM OF ACTION

UNIQUE MECHANISM OF ACTION


Floseal Matrix targets the beginning and end of the coagulation cascade: the **gelatin granules** assist with contact activation and platelet aggregation while the **thrombin** converts fibrinogen to fibrin to accelerate clot formation.¹³ The gelatin granules swell approximately 10-20% within 10 minutes upon contact with blood or other fluids creating a **tamponade effect**.¹⁰

PRODUCT APPLICATION¹⁰




APPLY

Identify the source of bleeding at the tissue surface and apply a small "mound" of **Floseal** Matrix to the target site. Do not inject or compress **Floseal** Matrix into blood vessels. Do not apply **Floseal** Matrix in the absence of blood flow.



APPROXIMATE

Immediately and gently apply approximate **Floseal** Matrix to the bleeding surface with a gauze sponge moistened with sterile (non-heparinized) saline for approximately 2 minutes.



IRRIGATE EXCESS

Once bleeding has ceased, excess **Floseal** Matrix not incorporated into the clot should always be removed by gentle irrigation and suctioned away.



LEAVE IN SITU

Do not disrupt the clot by physical manipulation. **Floseal** Matrix incorporated in the hemostatic clot should be left in situ. **Floseal** Matrix is biocompatible and resorbed within 6 to 8 weeks, consistent with normal wound healing.

FIRST AND ONLY¹⁰

Floseal Matrix is the **FIRST AND ONLY** active flowable hemostat to use recombinant DNA technology. The thrombin component is **Recothrom** Thrombin Topical (Recombinant) and has the following key benefits:



No human blood components



Structurally similar to human thrombin



Eliminate reliance on human blood donations

Selected Risk Information for FLOSEAL Matrix: Do not administer to patients with a history of or hypersensitivity to **Recothrom** Thrombin Topical [Recombinant], any component of **Recothrom** or hamster proteins.


PROVEN CLINICAL DATA^{10,15-17}

Study Design

Prospective, randomized controlled trials (RCT) have demonstrated that **Floseal** Matrix acts fast to achieve hemostasis, decreasing the potential for complications and reducing the procedure time. The objective of the pivotal studies was to evaluate the safety and effectiveness of **Floseal** Matrix compared to a **Gelfoam** Hemostat + Thrombin in controlling intraoperative bleeding. **Floseal** Matrix is effective across surgical specialties: 309 patients participated in the RCTs for procedures in cardiac, vascular, spinal or orthopedic surgeries.


Primary Endpoint: Cessation of Bleeding Within 10 Minutes

Hemostasis Within 10 Minutes First Lesion Only (Intent-to-Treat Patients)



Secondary Endpoint: Time To Hemostasis

Median Time to Hemostasis (Minutes)



Selected Risk Information for FLOSEAL Matrix: **Floseal** Matrix is not intended as a substitute for meticulous surgical technique and proper application of ligatures or other conventional procedures for hemostasis. **Floseal** Matrix is not intended as a prophylactic hemostatic agent.

PROVEN IN CARDIOVASCULAR SURGERY

FLOSEAL Matrix is effective in heparinized patients.^{10,15}



EFFECTIVE HEMOSTASIS IN HEPARINIZED CARDIAC SURGERY PATIENTS^{10,15}

Study Design


A prospective, randomized, multi-center trial comparing **Floseal Matrix** with a control hemostat (**Gelfoam Hemostat + Thrombin**) in stopping intraoperative bleeding during cardiovascular surgery was undertaken.

Four cardiac centers randomized 93 patients; 88 of 93 patients underwent extracorporeal cardiopulmonary bypass. Randomization occurred after bleeding requiring a topical hemostatic agent was identified. The primary endpoint was hemostasis for the first treated bleeding site within 10 minutes.

All patients were heparinized during the operation. The average heparin dose was $41,542 \pm 21,010$ U (mean $\pm$ SD) in the **Floseal Matrix** group and $38,289 \pm 15,506$ U (mean $\pm$ SD) in the control group.

The success rate for **Floseal Matrix** did not appear to be affected by whether or not the patient had received protamine sulfate administration. This was demonstrated by the success rate for **Floseal Matrix** before protamine sulfate administration which was similar to the success rate after protamine sulfate administration whereas the Control hemostat success rate was clearly lower before protamine sulfate reversal of heparin was administered.

There was no difference in the adverse event profile between the two groups.




REDUCTION IN BLOOD PRODUCT TRANSFUSIONS, SURGICAL REVISION RATES, OPERATIVE HEMOSTASIS TIME AND MINOR COMPLICATIONS^{10,14}

Study Design


In a prospective, randomized controlled trial for cardiovascular surgery, patients were assigned to a **Floseal Matrix** group (n=209) or the comparator group (n=206). The comparator group was treated with either **Surgicel Nu-Knit Absorbable Hemostat** (oxidized regenerated cellulose) or **Gelfoam** sponge (purified porcine skin gelatin).

A sub-analysis of patients who suffered from intraoperative bleeding (**Floseal Matrix** [n=110] and comparator [n=104]) revealed that **Floseal Matrix** is proven to reduce transfusion of blood products (p<0.001), reduce surgical revision rate for bleeding (p=0.04), reduce operative hemostasis time* (p<0.001), and reduce minor complications† (p=0.04).


Transfusion of Blood Products



Surgical Revision Rate for Bleeding




Operative Hemostasis Time*



*Operative time comprised between decannulation and closure of the sternum

Minor Complications†



†Defined as either renal failure, respiratory insufficiency, or inotropic support lasting more than 24 hours.

PROVEN IN
SPINE SURGERY

FLOSEAL Matrix is
2x faster*
at achieving hemostasis^{10,16}

*Compared to **Gelfoam** Hemostat + Thrombin.




The safety and effectiveness of **Floseal** Matrix in spinal surgery procedures have been proven.^{10,16}

A prospective, multi-center, randomized controlled trial comparing **Floseal** Matrix with **Gelfoam** Hemostat + Thrombin in stopping intraoperative bleeding during spinal surgery was undertaken. Four centers randomized 127 patients, 21 years of age or older. Randomization occurred after a bleeding lesion requiring a topical hemostatic agent was identified.^{10,16}

The primary endpoint was hemostasis for the first treated bleeding site. There was no difference in the adverse event profile between the two groups.^{10,16}

In this trial, the median time to hemostasis was 1.5 minutes (95% confidence interval [CI], 1, 1.5) for **Floseal** Matrix and 3 minutes (95% CI, 2, 4.5) for **Gelfoam** Hemostat + Thrombin.^{10,16}

GELFOAM Hemostat+Thrombin




Median time to hemostasis:
3.0 minutes^{10,16}



MORE EFFECTIVE AT EACH TIME POINT

Time to hemostasis in spinal surgery, first bleeding site^{10,16}




Floseal Matrix is More Effective than **Gelfoam** Hemostat + Thrombin at Each Time Point^{10,16}

At each time point, the percentage of success in the **Floseal** Matrix group was greater than the **Gelfoam** Hemostat + Thrombin group; this was also true for all bleeding sites ($p < 0.001$). The application was considered successful if the bleeding stopped within 10 minutes.^{10,16}

IMPROVED
VISUALIZATION

In spinal surgery procedures, excessive and/or uncontrolled bleeding imparts challenges related to adequate visualization.⁶



BEFORE
FLOSEAL MATRIX

Baxter Image on File



AFTER
FLOSEAL MATRIX

Baxter Image on File

Selected Risk Information for FLOSEAL Matrix: Do not inject or compress **Floseal** Matrix into blood vessels. Do not apply **Floseal** Matrix in the absence of active blood flow, e.g. while the vessel is clamped or bypassed, as extensive intravascular clotting and even death may result.

POTENTIAL FOR **SIZABLE COST SAVINGS** IN **CARDIAC SURGERY**


Study Design

Based on clinical outcomes from a prospective randomized clinical trial, Nasso et al¹⁴, a cost consequence framework was utilized to model the economic impact of comparator groups.¹⁸ Clinical outcomes were obtained and analyzed for a flowable hemostatic matrix (**Floseal** Matrix) vs passive, basic (non-flowable) topical hemostats (**Surgicel** and **Gelfoam** Hemostats).¹⁸

Costing analyses focused on the following outcomes: complications, blood transfusions, surgical revisions, and operating room (OR) time.¹⁸ The results suggest that if an active flowable hemostatic matrix (rather than a passive, basic hemostat) was utilized exclusively in 600 mixed cardiac surgeries annually, a hospital could improve outcomes by a reduction of 109 complications, 54 surgical revisions, 194 transfusions, and 242 hours of OR time.¹⁸

Study outcomes correspond to a **net annualized cost consequence savings of \$5.38 million when using FLOSEAL Matrix** versus passive hemostats (**Surgicel/Gelfoam** Hemostats) with complication avoidance as the largest contributor.¹⁸

Using **Floseal** Matrix Results in Significant Cost Savings¹⁸



*An avoidable cost consequence would result until the difference between **Floseal** Matrix and the **Surgicel/Gelfoam** reached \$8,960.36 per average case. At this price differential, there would no longer be an economic cost consequence between the comparative groups.

Clinical Outcomes **AVOIDED** by Using **Floseal** Matrix Instead of Passive Hemostats¹⁸

Outcome	Number Avoided
Major complications*	33 events
Minor complications	76 events
Surgical revisions for bleeding	54 events
Blood transfusions	194 events
Operating room time	242 hours

*Difference in rate of major complications between **Floseal** Matrix group and comparator group was not found to be statistically significant in the trial outcomes.

POTENTIAL FOR **SIZABLE COST SAVINGS** IN **SPINE SURGERY**

Study Design

Based on outcomes from Price et al⁷, a cost consequence framework was utilized to model the economic impact of **Floseal** Matrix vs. **Surgiflo** Matrix in spine surgery.^{7,19}

Compared to **Surgiflo** Matrix (base case of 182 major and 59 severe spine surgeries annually), **Floseal** Matrix required 3 fewer blood transfusions and reduced 53 hours of OR time, resulting in annual savings of \$151 per major and \$574 per severe spine surgery patient in addition to totally offsetting the product cost.^{7,19}

Study outcomes indicate that using **FLOSEAL Matrix** instead of **SURGIFLO** Matrix in spine surgery potentially **leads to sizable cost savings for hospitals**. Monte-Carlo simulations demonstrated **annual savings in major and severe spine surgery were \$27,412 and \$33,839, respectively**.¹⁹

Economic Outcomes for **Floseal** Matrix versus **Surgiflo** Matrix
Base Case: Medium-Size Hospital*^{7,19}

MAJOR SPINE SURGERY	
Medium-Size Hospital, Mean Number of Surgeries (Range)*	182 Major Spine Surgeries (101-300)
Number of OR Hours Saved	27
Number of Transfusions Avoided	3 events
Net Cost Savings†	\$27,412 (\$151 Per Surgery)
SEVERE SPINE SURGERY	
Medium-Size Hospital, Mean Number of Surgeries (Range)*	59 Severe Spine Surgeries (39-100)
Number of OR Hours Saved	26
Net Cost Savings†	\$33,389 (\$574 Per Surgery)

Outcomes associated with **FLOSEAL** Matrix correspond to an **ANNUAL SAVINGS OF**

\$27,412


IN MAJOR SPINE SURGERIES AND

\$33,839

IN SEVERE SPINE SURGERIES

*In an analysis of all spinal surgeries performed in US hospitals in 2012 HCUP sample of n=2,071 hospitals, a medium-volume hospital represents the 51st-82nd percentile. (Low-volume hospitals were below 50th percentile while high-volume hospitals represented 83rd percentile and higher.) The mean number of surgeries per hospital in each group was then calculated for use in the cost model.

†Based on 2016 US Dollar.



Ordering Information

PRODUCT	QTY	PRODUCT NUMBER
FLOSEAL Matrix 5mL	Case of 6	ADS202105
FLOSEAL Matrix 10mL	Case of 6	ADS202110
Endoscopic Applicator (41 cm)	Pack of 6	0600225
13 cm Malleable Applicator Tip	Pack of 6	ADS201816
Disposable Curved Applicator Tip	Pack of 10	ADS201815

FLOSEAL Hemostatic Matrix Indications and Important Risk Information

FLOSEAL Hemostatic Matrix Indication

Floseal Matrix is indicated in surgical procedures (other than ophthalmic) as an adjunct to hemostasis when control of bleeding by ligature or conventional procedures is ineffective or impractical.

Important Risk Information

Do not inject intravascularly.

Do not inject or compress Floseal Matrix into blood vessels.

Do not apply Floseal Matrix in the absence of active blood flow, e.g., while the vessel is clamped or bypassed, as extensive intravascular clotting and even death may result.

Do not administer to patients with a history of hypersensitivity to Recothrom Thrombin topical (Recombinant), any components of Recothrom, or hamster proteins.

Do not use Floseal Matrix in patients with known allergies to materials of bovine origin.

Do not use Floseal Matrix in the closure of skin incisions because it may interfere with the healing of the skin edges due to mechanical interposition of gelatin.

Floseal Matrix is not intended as a substitute for meticulous surgical technique and the proper application of ligatures or other conventional procedures for hemostasis. Floseal Matrix is not intended to be used as a prophylactic hemostatic agent.

As with other hemostatic agents, do not apply Floseal Matrix to sites where there is negative peripheral venous pressure, as material may be drawn into the vascular system potentially resulting in life-threatening thromboembolic events.

Excess Floseal Matrix (material not incorporated in the hemostatic clot) should always be removed by gentle irrigation and suctioned out of the wound.

The particles of Floseal Matrix swell approximately 10-20% (upon contact with blood or other fluids) and surgeons should consider its potential effect on the surrounding anatomic areas. Maximum swell volume is achieved within about 10 minutes.

Floseal Matrix should not be used in the presence of infection. Floseal Matrix should be used with caution in contaminated areas of the body.

The safety and effectiveness of Floseal Matrix for use in ophthalmic procedures has not been established.

Floseal Matrix should not be used for controlling post-partum bleeding or menorrhagia.

The safety and effectiveness of Floseal Matrix has not been established in children under 2 years of age and pregnant women.

It is not known whether Floseal Matrix can cause fetal harm when administered to a pregnant woman or can affect reproduction capacity. Floseal Matrix should be administered to a pregnant woman only if clearly needed.

Do not use air to remove residual Floseal Matrix from Applicator tip. The Applicator tips should not be cut, as tissue injury from sharp edges may result.

Hypersensitivity reactions, including anaphylaxis, may occur. Recothrom thrombin is produced in a genetically modified Chinese Hamster Ovary (CHO) cell line and may contain hamster or snake proteins.

For single use only. Do not re-sterilize.

Floseal Matrix should not be applied before the application site is cleaned to remove any antiseptics that may contain alcohol, iodine, or heavy metal ions.

When placed into cavities or closed tissue spaces, gentle approximation is advised. Do not compress.

As with other hemostatic agents, do not aspirate Floseal Matrix into extracorporeal cardiopulmonary bypass circuits or autologous blood salvage circuits.

Do not use Floseal Matrix on bone surfaces where adhesives, such as methylmethacrylate or other acrylic adhesives, will be required to attach a prosthetic device.

Floseal Matrix should not be used for the primary treatment of coagulation disorders.

The safety and effectiveness of the combined use of Floseal Matrix with antibiotic solutions or powders has not been established.

The safety and effectiveness for use in neurosurgical and urological procedures has not been established through randomized clinical studies.

In urological procedures, Floseal Matrix should not be left in the renal pelvis or ureters to eliminate the potential foci for calculus formation.

Antibody formation to Recothrom thrombin occurred in <1% of patients. None of the antibodies detected neutralized native human thrombin.

Thrombin must be added to the Gelatin Matrix prior to use.

Rx Only. For safe and proper use of this device, refer to the full Instructions for Use.

References

- Iannitti DA, Kim C, Ito D, et al. Impact of an active hemostatic product treatment approach on bleeding-related complications and hospital costs among inpatient surgeries in the United States. J Med Econ. 2021;24(11):514-523.
- Stokes ME, Ye X, Shah M, et al. Impact of bleeding-related complications and/or blood product transfusions on hospital costs in inpatient surgical patients. BMC Health Serv Res. 2011; 11:135.
- Preliminary Results, Premier Data Summary, Report_2018Q1_2019Q2 (v1.0). Baxter [Data on File]. March 9, 2020. Study number 9.4.
- Centers for Disease Control and Prevention. Chronic Liver Disease and Cirrhosis. <https://www.cdc.gov/nchs/fastats/liver-disease.htm>.
- Stuntz M, Bernstein B. Recent trends in the prevalence of low-dose aspirin use for primary and secondary prevention of cardiovascular disease in the United States, 2012-2015. Prev Med Rep. 2017;5:183-186.
- Ramirez MG, Deutsch H, Khanna N, et al. Floseal only versus in combination in spine surgery: a comparative, retrospective hospital database evaluation of clinical and healthcare resource outcomes. Hosp Pract. 2018;46(4):189-196.
- Price JS, Tackett S, Patel V. Observational evaluation of outcomes and resource utilization from hemostatic matrices in spine surgery. J Med Econ. 2015;18(10):777-786.
- Lewis KM, Li Q, Jones DS, et al. Development and validation of an intraoperative bleeding severity scale for use in clinical studies of hemostatic agents. Surgery. 2017;161:771-781.
- Efficacy and Safety of Active Hemostatic Agents at Specific Bleeding Grades as Determined by the Validated Intraoperative Bleeding Scale (VIBe Scale) in Heparinized Porcine Bleeding Models. Baxter [Data on File]. Study number BWQ026-IS17/BWQ027-IS17.
- Floseal Hemostatic Matrix Instructions for use.
- Lewis KM, Atlee HD, Mannone AJ, et al. Comparison of two gelatin and thrombin combination hemostats in a porcine liver abrasion model. J Invest Surg. 2013;26:141-148.
- Floseal Matrix with Recothrom Thrombin topical [Recombinant] Time and Motion Study. Baxter [Data on File]. Study number BXU569412. March 25, 2022.
- Oz MC, Rondinone JF, Shargill NS. Floseal Matrix: new generation topical hemostatic sealant. J Card Surg. 2003;18:486-493.
- Nasso G, Piancone F, Bonifazi R, et al. Prospective, randomized clinical trial of the Floseal Matrix sealant in cardiac surgery. Ann Thorac Surg. 2009;88(5):1520-1526.
- Oz MC, Cosgrove DM, Badduke BR, et al. Controlled clinical trial of a novel hemostatic agent in cardiac surgery. Ann Thorac Surg. 2000;69(5):1376-82.
- Renkens KL, Payner TD, Leipzig TJ, et al. A multicenter, prospective, randomized trial evaluating a new hemostatic agent for spinal surgery. Spine. 2001;26(15):1645-50.
- Weaver FA, Hood DB, Zatina M, et al. Gelatin-thrombin-based hemostatic sealant for intraoperative bleeding in vascular surgery. Ann Vasc Surg. 2002;16(3):286-93.
- Tackett SM, Sugarman R, Kreuwel HT, et al. Hospital economic impact from hemostatic matrix usage in cardiac surgery. J Med Econ. 2014;17(9):670-676.
- Makhija D, Rock M, Ikeme S, et al. Cost consequence analysis of two different active flowable hemostatic matrices in spine surgery patients. J Med Econ. 2017;12:1-8.

advancedsurgery.baxter.com

Customer Service: 1-888-229-0001

Baxter, Floseal and Recothrom are trademarks of Baxter International Inc. or its subsidiaries.

US-AS19-250002 (v1.0) 04/2025

Baxter